

Loupe selector

Step 1 – Select your magnification

Discipline	Magnification
Cardiology	2.5x – 3.0x
ENT	2.5x – 3.0x
Neurosurgery	2.5x – 3.0x
Ob. Gyn.	2.5x – 3.0x
Ophthalmology	2.5x – 3.0x
Orthopaedic	2.5x – 3.0x
Plastics	2.5x – 3.0x
Genito-Urinary	2.5x – 3.0x
Surgical assistants	2.5x – 3.5x
Dental – General	2.5x
Dental – Endo, crown, bridge	3.5x – 4.5x
Dental – Lab work, technician	4.5x – 5.0x

Step 2 – Select your working distance

Height	Sitting	Standing
< 5' 7" 150cm - 170cm	34cm	34cm
5' 7" – 5' 10" 170cm – 178cm	34cm	42cm
5' 10" – 6' 4" 178cm - 193cm	42cm	46cm
> 6' 4" 193cm	46cm	50cm

Step 3 – Select your frame

Distance between temples	Frame selection
130mm	Small frame
135mm	Regular frame
135mm	LT (Long Temple) frame 150mm arm length